



CE

**NOVO !****RANDOX** REAGENTS  
Cardiology**TxB** Cardio™**NOVINE U DIJAGNOSTICI  
U OBLASTI KARDIOLOGIJE****SUPERLAB** Vam u saradnji sa Randox kompanijom predstavlja novi proizvod, **TxB reagens** !**TxB Cardio- reagens za utvrđivanje rezistencije na aspirin**

**Aspirin** je antiagregacioni lek koji se široko primenjuje u primarnoj i sekundarnoj terapiji kardiovaskularnih bolesti. Pacijenti koji na adekvatan način ne reaguju na terapiju aspirinom imaju mnogo veći rizik od smrti zbog neželjenih kardiovaskularnih događaja od ljudi sa pozitivnim odgovorom na terapiju. Identifikacija ovih pacijenata je značajno unapređena upotrebom **TxB Cardio Randox reagensa**. Rezultati dobijeni ovim testom omogućavaju pravovremenu reakciju kliničara kada su u pitanju pacijenti iz visoko rizične grupe.

Randox TxB Cardio Assay  
Details**Methodology**

Latex-enhanced immunoturbidimetric

**Sample Type**

Urine

**Reagent Format**

Liquid ready-to-use

**Expiry**Stable to expiry date when stored at  
+2 to +8°C**Assay Range**

400-6000 pg/ml

**QC Material**

Three levels of controls available

**Calibrator Information**

TxB Cardio Calibrator Series available

**Cat. No.:**TBX2759 R1 1 x 9ml  
R2 1 x 4,7ml**Da li ste znali da 25-30% pacijenata koji primaju niske doze aspirina u terapiji pati od rezistencije na aspirin?**

- Rezultati analize pacijenata 20 studija pokazali su da 28% u grupi od 2930 ispitivanih je rezistentno na aspirin.
- Aspirin rezistentne osobe pokazuju 2x veći rizik od kardiovaskularnih događaja i 3.5x veći rizik od smrti.
- **TxB Cardio** je latex imunoturbidimetrijski test sa instrument specifičnim aplikacijama, tako da ne zahteva posebnu opremu za izvođenje.

**RANDOX**Randox Laboratories Limited, 55 Diamond Road, Crumlin, County Antrim, BT29 4QY, United Kingdom  
T +44 (0) 28 9442 2413 F +44 (0) 28 9445 2912 E marketing@randox.com I www.randox.com

Information correct at time of print. Randox Laboratories Limited is a subsidiary of Randox Holdings Limited, a company registered within Northern Ireland with company number NI 614690. VAT Registered Number: GB 151 6827 08. Product availability may vary from country to country. Please contact your local Randox representative for information. Products may be for Research Use Only and not for use in diagnostic procedures in the USA.



CE

L276 (09/11)

**SUPERLAB**®  
Your Lab - Our PassionM. Milankovića 25, Tel./Fax: 011 22 22 222, lokal 275  
E-mail: bkrivokuca@super-lab.com www.super-lab.com